

Fabric Thickness Tester II™(Digital)

INTRODUCTION

Determination of thickness of fabric samples in laboratory is usually carried out with the help of a precision thickness gauge. In this equipment, the fabric whose thickness is to be determined is kept on a flat anvil and a circular pressure foot is pressed on to it from the top under a standard fixed load. The Digital Display gives the Thickness in mm.

EQUIPMENT

Paramount **Fabric Thickness Tester II™** main body is completely made of mild steel. Digital Gauge for measuring the thickness is mounted on top. Specially designed stainless steel anvil is provided to keep the sample to be tested. Two different size of pressure foot are provided for testing different types of fabrics. A standard load is mounted on the dial to give standard pressure on the sample being tested. The equipment is meticulously finished in blue metallic paint.

SALIENT FEATURES

- Specialized equipment to determine the thickness of fabrics digitally.
- Latest model with modern aesthetics.
- Maximum capacity 10mm & accuracy 0.01mm.
- Portable & handy to carry anywhere.
- Smooth precision-engineered components for accurate results.
- Analog dial for measuring thickness in millimeters.
- Supplied with two pressure feet of different size for different type of fabrics.
- A standard Dead weight is provided with the instruments for accurate result.
- Supplied with calibration & inspection certificates.

Fabric Thickness Tester™ is supplied complete with below accessories.

Main Unit (With 10 mm Pressure Foot Fitted)	: 01 No.
Pressure Foot (Dia. 25 mm)	: 01 No.
Pressure Foot (Dia. 10 mm)	: 01 No.
This User's Manual	: 01 No.
Inspection & Conformance Certificate	: 01 No.
Calibration Certificate (Traceable to NPL)	: 01 No.



TECHNICAL DETAILS

Range of measurement	: 0 - 10 mm
Least count	: 0.01 mm
Diameter of anvil	: 60 mm
Diameter of pressure foot (Interchangeable)	: 10 mm & 25 mm
Load on pressure foot	: 78.5 g (100g/sq. cm)
Throat depth	: 22 mm
Overall Dimensions of the Instrument	: 209.5 (H) x 152 (D) mm 8.25" (H) x 6" (D) Inch
Net Weight of the Instrument	: 3.4 Kg. (7.5 lbs)

RELATED STANDARDS

- ASTM D 1777 – 96
- ASTM D 5729 – 97
- ISO 3616/5084/9073



Corporate Member



An ISO 14001 : 2004 Group
 We Care for our Environment !



Member : Bureau of Indian Standards