

High Precision Digital Balance

INTRODUCTION

Paramount *high precision digital balance* is very useful equipment to determine count of yarn & GSM of fabrics with high accuracy. This high precision balance weighs up to four decimal numbers, so it is mostly used in textile laboratories, College & Universities to weigh small sample of yarn & fabrics.

EQUIPMENT

Paramount *high precision digital balance* is very delicate weighing balance with high accuracy up to four decimal numbers. Specially designed acrylic cover is provided for minimizing the external interferences while weighing. A.C adaptor is supplied with the equipment to run in 240V A.C supply.

SALIENT FEATURES

- Specially designed high precision digital balance used to determine GSM of any fabrics & Count of Yarns with high accuracy.
- High precision digital balance high accuracy up to four decimal numbers.
- Maximum capacity 210gms.
- With LCD display.
- Sprit level for leveling of the equipment.
- Latest model with Auto calibration feature.
- Supplied with inspection & calibration certificates




High Precision Digital Balance is supplied complete with below accessories.

Main unit .	: 01 No.
Acrylic Cover	: 01 No.
Adaptor	: 01 No.
This User's Manual	: 01 No.
Inspection & Conformance Certificates	: 01 No.
Calibration Certificates	: 01 No.

TECHNICAL DETAILS

Capacity	:210 gms
Sensitivity	:0.0001gms
Pan	:stainless steel
Cover (Air)	:acrylic
Sprit Level	:For leveling
Input	:240V
Dimension	: 260(L)x, 165(W)x, 100(H)mm 10(L)x, 6.5(W)x, 3.9(H)inch
_Weight	: 1.2Kg (2.6 4 lbs)



AATCC
CORPORATE MEMBER



Patron Member



Member Bureau of Indian Standards