

*swatch*CUTTER 500 M™(Eco Model)

INTRODUCTION

Paramount *swatch*CUTTER 500 M™(Eco Model) is an economical swatch cutting machine for preparing swatches with pinking edges in multiple layers up to a depth of 40 mm (1.5") and width of 500 mm (19.5).

EQUIPMENT DETAILS

Paramount *swatch*CUTTER 500 M™ (Eco Model) is completely made out of mild steel. The cutting action is accomplished by specially imported Italian blade with the help of Motorized System. The cutting table itself moves upward to cut the sample while the blade is fixed. A special grade cutting pad is provided to help in cutting swatches & for the long life of the blade. One scale is mounted in one side of the machine for measuring the sample length to cut. The equipment is mounted on 4 Custer wheels with braking system to ease in transportation.

SALIENT FEATURES

- An economical motorized swatch cutting machine for making samples as per size & design.
- Sample of maximum 40cm width & any length can be prepared.
- Economically design with digital control.
- Smooth precision-engineered components for exceptional performance.
- Pinking blade tip gives excellent look of swatches prepared.
- Caster wheel at base with stopper gives flexibility in transport.
- A Calibrated scale attached to prepare sample as size required.
- Supplied with inspection certificates.

***swatch*CUTTER 500 M™(Eco Model) is supplied complete with below accessories.**

Main Unit	: 01 No.
Special Cutting Mat / Board.	: 01 No.
User's Manual	: 01 No.
Inspection & Conformance Certificate.	: 01 No.

TECHNICAL DETAILS

Electrical supply	: 230V, Single Phase AC, 50Hz
Motor	: 1-φ AC, 230 volt 50Hz. 0.5HP, 1400 rpm.
Cutting Width	: 500mm.
Cutting Height	: 40mm.
Dimension of the cutting pad	: 520mm x 520mm
Graduated scale	: 50 cm left & right both side from 'O' Position / cut position on the both edge. L/c: 1mm.
Dimension of the instrument	: 952mm (L) x 660mm (D) x 1140mm (H). 37.5 Inch (L) x 26 Inch (D) x 43.5 inch (H)
Net weight of the instrument	: 135 kg (297 lbs)



**AATCC
CORPORATE MEMBER**



Patron Member



**Member Bureau of
Indian Standards**