

swatchCUTTER 500™
(Fabric Sampling (Pinking) Machine)

INTRODUCTION

PARAMOUNT swatchCUTTER 500™ is an accurate & convenient Equipment for preparing swatches with pinking edges in multiple layers to a depth of 20 mm (0.8") and width of 500 mm (19.5)

EQUIPMENT

PARAMOUNT swatchCUTTER 500™ is a user friendly versatile Swatch Cutting Machine for making swatches upto 500mm (19.5") quickly & conveniently The specially designed Equipment can cut swatches in multiple layers up to a depth of 20mm. It can make samples out of any type of fabric (woven/knitted or non woven) with the special imported Italian scissor. The equipment is mounted on a strong & sturdy table for ease in operation. The sample sizes can conveniently determined by the graduated screen fixed on the fabric guide bar. The cutting area is fully protected by a transparent acrylic sheet for safety, which also helps the operator to view the cutting operation. The specially designed cutting board is provided to support the sample to cut & long life of the blade. A specially designed long handle is attached with the equipment for smooth cutting with minimum operator fatigue. The equipment is meticulously finished in blue & grey paint to give a gorgeous look.

SALIENT FEATURES

- Specially designed equipment to prepare swatches of different size & design.
- Sample up to 500mm width & any length can be cut.
- Depth of the sample up to 20mm can easily be cut.
- Wooden stand to keep the instrument stable.
- Smooth precision-engineered components for exceptional performance.
- The pinking blade tip gives excellent look of swatches prepared.
- A Calibrated scale attached to prepare sample of size required.
- Supplied with inspection certificates.

swatchCUTTER 500™ (Normal) is supplied complete with below accessories.

Main Unit	: 01 No.
Sturdy Table (Fitted or Dismantled Position)	: 01 No.
Special Cutting Mat/Board (Made in Italy)	: 01 No.
Double End Spanner - (13-13)	: 01 No.
Double End Spanner - (10-10)	: 01 No.
This User's Manual	: 01 No.
Calibration Certificate (Traceable to NPL)	: 01 No.
Inspection & Conformance Certificate	: 01 No.



TECHNICAL DETAILS

Cutting Width	: 500 mm
Cutting Height	: 20 mm
Dimension of the cutting pad	: 500 x 500 mm (±2 mm)
Graduated Scale	: 55.5 Cm. on left side & 38.5 Cm. on right
	side; From the "0" position (cut position)
	L/C : 1 mm
Dimension of the Instrument (Including Table)	: 1015 (L) x 775 (W) x 1370 (H) mm 28" (W) x 24" (D) x 52" (H) Inch.
Dimension of the Instrument (Excluding Table)	: 950 (L) x 790 (W) x 600 (H) mm 38" (W) x 32" (D) x 29" (H) Inch.
Net Weight of the Instrument	: 45 Kg.
Net Weight of the Table	: 42 Kg



AATCC
CORPORATE MEMBER



Patron Member



Member Bureau of Indian Standards