

# TDS Meter

**INTRODUCTION**

Content of Total Dissolved Salt is an essential phenomenon for performing chemical test in the laboratory. **Paramount TDS Meter** is a handy instrument to instantly check TDS value of any liquid in PPM.

**EQUIPMENT**

**Paramount TDS Meter** is made out of high quality plastic material. It is handy equipment having a sensing probe to sense the TDS in water. The Range is 0 to 1999 PPM. The equipment is specially sealed to protect from water to get damaged.

**SALIENT FEATURES**

- Used to evaluate the TDS (Total Dissolved Salts) in PPM of water.
- Latest design with modern aesthetics.
- Portable in Size
- Digital display of TDS value.
- Plastic Body for lightweight & flexibility.

**TDS Meter is supplied complete with below accessories.**

|                                      |          |
|--------------------------------------|----------|
| Main Unit                            | : 01 No  |
| This User's Manual                   | : 01 No. |
| Inspection & Conformance Certificate | : 01 No. |

**TECHNICAL DETAILS**

|               |                             |
|---------------|-----------------------------|
| Range         | : 0 -1999 PPM               |
| Least Count   | : 1                         |
| <b>Weight</b> | <b>: 88 Gm (0.1936 lbs)</b> |

**Dimensions:**

|        |          |
|--------|----------|
| Length | Diameter |
| 175 mm | 40 mm    |
| 7 inch | 1.5 inch |



**AATCC**  
CORPORATE MEMBER



**Patron Member**



**Member Bureau of Indian Standards**