

textile LabSHAKER

INTRODUCTION

Paramount *textile LabSHAKER* is very useful equipment for Textile College, Universities & laboratories for the preparation of dye solution & chemical solution. The equipment helps in shaking or agitating the solution/liquor at predetermined speed for desired span of time.

EQUIPMENT

Paramount *textile LabSHAKER* main body is made of mild steel with a top-rotating unit, which can accommodate 16 Conical flask or beakers. The Specially designed top rotating unit having twenty-five rubber diaphragms attached with the help of screws for holding the beakers without slippage. The speed of the rotating unit can be varied with speed regulator knob. The equipment is meticulously finished in cream colour paint.

SALIENT FEATURES

- Specially designed equipment to be used in laboratories for continuous Shaking of liquor / liquid.
- With speed regulation system to change the level of shaking.
- Up to 16 Conical flask or beakers can be kept together for shaking.
- Specially designed rubber beaker holders give flexibility to use glass Conical flask or beakers without damage.
- Ergonomically designed with modern aesthetics.
- Smooth precision-engineered components for exceptional performance.
- Supplied with all accessories.
- Supplied with Inspection & Conformance Certificate.

textile LabSHAKER is supplied complete with below accessories

Main Unit	: 01 No.
Dust Cover	: 01 No.
500 ml Conical Flask	: 04 No.
Detailed User's Manual	: 01 No.
Calibration Certificate (Traceable to NPL)	: 01 No.
Inspection & Conformance Certificate	: 01 No.

TECHNICAL DETAILS

Dimension:	: 508(L) x, 508(W) x, 375(H) mm
	: 20(L) x, 20(W) x, 14.75(H) inch
Weight	: 16Kg (35.02lbs)



AATCC
CORPORATE MEMBER



Patron Member



Member Bureau of
Indian Standards